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
May 2016

Information Regarding
LIFEPAK® 15 MONITOR/DEFIBRILLATOR WITH End-tidal CO₂ OPTION

Products: ~ This information only applies to LIFEPAK® 15 monitor/defibrillators with the End-tidal CO₂ (EtCO₂) option configured to the kPa or % mode. If you use EtCO₂ in mmHg this is not applicable.

Dear Valued Customer,

Physio-Control, Inc. LIFEPAK 15 monitor/defibrillators with EtCO₂ capability are distributed with the EtCO₂ setting preset to display in mmHg. Customers may choose to reconfigure the device to display the EtCO₂ in alternative modes of kPa or % settings. We have become aware that when using EtCO₂ in the kPa or % setting and in a situation where the reading is above 9.9 kPa, the display of the LIFEPAK 15 respiratory rate may partially obscure a portion of the leading digit of the EtCO₂ value. See examples on page 2.

Although part of the first digit of the EtCO₂ value is obscured by the respiratory rate, the capnography waveform scale automatically adapts to the EtCO₂ value. Observing the capnography waveform and the range displayed will help the user confirm the patient's EtCO₂ value.

Clinical interpretation of a low EtCO₂ value when the value is actually high could potentially lead to incorrect treatment. There have been no reports of incorrect patient treatment or harm associated with this issue.

The LIFEPAK 15 Operating Instructions for monitoring EtCO₂ include displaying the capnography waveform and instructions for setting high and low EtCO₂ alarms. Along with patient assessment, the user should display the EtCO₂ waveform and activate the alarms to monitor EtCO₂ values.

Actions:

1. Please forward this information to all of your sites, trainers and users that may have a LIFEPAK 15 monitor/defibrillator.
2. If you use your LIFEPAK 15 monitor/defibrillator to display EtCO₂ values you should also display the EtCO₂ waveform which includes the range along the right side of the display.
3. If you use either kPa or % settings, activate EtCO₂ alarms to identify EtCO₂ values outside the normal range.
4. If your protocols allow this option, consider changing the EtCO₂ to measure in the mmHg setting.

Physio-Control will post information on our website when updated software is available. If you want the update follow the instructions to update your device.

Should you have any questions about this subject, please contact us at 1-800-442-1142, option 7, 6:00 A.M. to 4:00 P.M. (Pacific), Monday - Friday.


Sincerely,

Rod J. Rylands
Vice President, Quality
PHYSIO-CONTROL, INC.

Physio-Control Nederland | Lifesaving starts here.™

LIFEPAK 15 showing EtCO₂ in % (percentage). The left side of the display shows the readings first character partially covered. The right side of the display also shows the range of readings and is showing that the level is well above 10%.

LIFEPAK 15 showing EtCO₂ in kPa (kilopascals). The left side of the display shows the readings first character partially covered. The right side of the display also shows the range of readings and is showing that the level is well above 10 kPa.