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BIO-RAD

Liquichek™ Diabetes Control Levels 1, 2 and 3

REF	171	Level 1	6 x 1 mL
	172	Level 2	6 x 1 mL
	173	Level 3	6 x 1 mL
	172X	Trilevel MiniPak	3 x 1 mL

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SUMMARY AND PRINCIPLE

The use of quality control materials is indicated as an objective assessment of the precision of methods and techniques in use and is an integral part of good laboratory practices. Three levels of control are available to allow performance monitoring within the clinical range. For customers in Germany, Quality control materials are required for assessment of laboratory performance as described in the "Guideline for Quality Assurance of Medical Laboratory Examinations following the German Medical Association" (RiL-BÄK regulation).

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STORAGE AND STABILITY

This product will be stable until the expiration date when stored unopened at -10 to -70°C. The control may be stored unopened at 2 to 8°C for 6 months, but should not be used past the expiration date (note the date refrigerated temperature storage begins). Once the control is opened, all analytes will be stable for 14 days when stored tightly capped at 2 to 8°C. Once thawed, do not refreeze the control; discard the remaining material.

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LIMITATIONS

1. This product should not be used past the expiration date.
2. If there is evidence of microbial contamination or excessive turbidity in the product, discard the vial.
3. This product is not intended for use as a standard.
4. ~~Wichtige Änderungen sind farblich hervorgehoben~~

ASSIGNMENT OF VALUES

The mean values printed in this insert were derived from replicate analyses and are specific for this lot of product. The tests listed were performed by the manufacturer and/or independent laboratories using manufacturer supported reagents and a representative sampling of this lot of control. Individual laboratory means should fall within the corresponding acceptable range, however, laboratory means may vary from the listed values during the life of this control. Variations over time and between laboratories may be caused by differences in laboratory technique, instrumentation and reagents, or by manufacturer test method modifications. It is recommended that each laboratory establish its own means and acceptable ranges and use those provided only as guides.

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[i] Significant changes are highlighted

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ENGLISH

Biological source material. Treat as potentially infectious.

Each whole blood donor unit used to manufacture this control was tested by FDA accepted methods and found non-reactive for Hepatitis B Surface Antigen (HBsAg), antibody to Hepatitis C (HCV) and antibody to HIV-1/HIV-2. This product may also contain other human source material for which there are no approved tests. In accordance with good laboratory practice, all human source material should be considered potentially infectious and handled with the same precautions used with patient specimens.

Safety Data Sheet (SDS) available for professional users on www.bio-rad.com.

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ENGLISH

- (1) NGSP units = %; IFCC unit = mmol/mol
 ▲ Data is not available at this time. Please Inquire.
 § The data required to establish the means and acceptable ranges for this assay were not obtained due to limited assignment participation. If your facility is interested in participating in the Value Assignment Program for this assay, please contact your local Bio-Rad office.
 ◆ INTERNATIONAL USE ONLY - The following section contains data for methods that are not available for diagnostic use in the United States.

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- INTERNATIONAL USE ONLY -

The following section contains data for methods that are not available for diagnostic use in the United States. ❖

FOR INTERNATIONAL USE ONLY

INTERNATIONAL USE ONLY -

The following section contains data for methods that are not available for diagnostic use in the United States.

METHOD

Methode // Méthode // Metodo // Método // Metodo // Metod // Metode // Yöntem // 測定法

	Units ⁽¹⁾	Level 1 - 38531		Level 2 - 38532		Level 3 - 38533		SI ⁽²⁾	Level 1 - 38531		Level 2 - 38532		Level 3 - 38533	
		Mean	Range	Mean	Range	Mean	Range		Mean	Range	Mean	Range	Mean	Range
HEMOGLOBIN A1C														
ARKRAY ADAMS A1C HA-8180 (Ver. 4.00) / HA-8170 / HA-8180 / HA-8180T / HA-8180V / HA-8181 / HYBRID AH-8280 / HYBRID AH-8290	%	5.00	4.50 - 5.50	9.30	8.40 - 10.2	13.9	13.0 - 14.8	mmol/mol	31.2	28.0 - 34.3	78.1	70.9 - 89.0	128	119 - 137
Bio-Rad VARIANT II HbA1c Dual Program 270-2500	%	▲		▲		▲		mmol/mol	▲		▲		▲	
Bio-Rad VARIANT II TURBO HbA1c Kit 2.0 270-2455 EX	%	▲		▲		▲		mmol/mol	▲		▲		▲	
NRIFON KONDEN Celltac Chemi HbA1c	%	5.20	4.50 - 5.90	9.80	8.40 - 11.2		12.1 - >13.0	mmol/mol	33.3	28.8 - 37.4	82.5	71.7 - 95.6		109 - >119
SEKISUI NORUDIA N HbA1c	%	4.41	3.53 - 5.29	8.59	6.87 - 10.3	13.2	10.6 - 15.8	mmol/mol	24.7	19.8 - 31.2	61.4	56.3 - 84.5	121	96.4 - 145
TOSOH GX Automated HPLC Analyzer	%	5.02	4.62 - 5.37	9.19	8.45 - 9.83	13.4	12.4 - 14.4	mmol/mol	31.4	28.8 - 34.3	78.9	70.8 - 83.1	123	113 - 133
TOSOH G11 Automated HPLC Analyzer (STD Mode)	%	4.83	4.44 - 5.17	8.90	8.19 - 9.52	13.0	12.0 - 13.9	mmol/mol	29.2	26.9 - 31.6	72.5	67.9 - 79.7	119	109 - 128
TOSOH HLC-723G8 / G9 / G11	%	5.20	4.70 - 5.70	8.70	8.70 - 10.7	14.5	13.5 - 15.9	mmol/mol	33.3	30.0 - 36.6	82.5	74.3 - 90.8	135	121 - 148

FOR INTERNATIONAL USE

NOTES

Bio-Rad Laboratories comprehensive line of quality controls and QC data management solutions.

Autoimmune Controls
Blood Gas Controls
Cardiac Assessment Controls
Chemistry Controls
Coagulation Controls
Congenital/Pediatric Disease Controls

Diabetes/Hemoglobin Controls
Hematology Controls
Hepatitis & HIV Tests Controls
Immunology Controls
Immunology/Protein Controls
Molecular Controls
Sexually Transmitted Disease Controls

Specialty Infectious Disease Controls
Therapeutic Drug Monitoring Controls
Toxicology: Drugs-of-Abuse Controls
Toxicology: Specialty Controls
Urinalysis Controls
QC Data Management Solutions
External Quality Assurance Services (EQAS)

BIO-RAD

USA: 1-800-824-6722 | www.bio-rad.com/qualitycontrol

MANUFACTURERS LISTED

Herstellerliste // Liste des fabricants // Elenco Produttori // Lista de fabricantes // Fabricantes enunciados // Lista över tillverkare // Liste over producenter // Üretici Listesi // 製造元一覧

Abbott Laboratories, Abbott Park, Illinois
 Arkray, Kyoto, Japan
 Beckman Coulter, Inc., Brea, California
 Beckman Coulter, Inc., Irving, Texas
 Bio-Rad Laboratories, Inc., Hercules, California
 Bio-Rad Laboratories, GmbH, Munich, Germany
 HemoCue, Ängelholm, Sweden
 Infopia, Gyeonggi-do, South Korea
 Nihon Kohden Corporation, Tokyo, Japan
 Ortho-Clinical Diagnostics, A Johnson and Johnson Company, Rochester, New York

Pointe Scientific, Inc., Canton, Michigan
 PTS Diagnostics, Indianapolis, Indiana
 Roche Diagnostics, GmbH, Penzberg, Germany
 Roche Diagnostics, Indianapolis, Indiana
 Sekisui Chemical Co. LTD., Tokyo, Japan
 Siemens Healthcare Diagnostics Inc., Newark, Delaware
 Siemens Healthcare Diagnostics Inc., Tarrytown, New York
 TOSOH Bioscience Inc., South San Francisco, California
 TOSOH Bioscience Inc., Tokyo, Japan
 Trinity Biotech, Kansas City, Missouri

UNITED STATES, Bio-Rad Laboratories
 9500 Jeronimo Road, Irvine, CA 92618

FRANCE, Bio-Rad
 3 boulevard Raymond Poincaré
 92430 Marnes-la-Coquette
 Phone: (33) 1-4793-6000 / Fax: (33) 1-4741-9133

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
UNITED STATES, Bio-Rad Laboratories
9500 Jeronimo Road, Irvine, CA 92618

FRANCE, Bio-Rad
3 boulevard Raymond Poincaré
92430 Marnes-la-Coquette
Phone: (33) 1-4795-6000 / Fax: (33) 1-4741-9133

**Bio-Rad
Laboratories**

Irvine • California • 92618 ☎ (800) 854-6737 • Telefax (949) 598-1550

BIO-RAD

Liquichek™ Diabetes Control

Levels 1, 2 and 3

REF	171	Level 1	6 x 1 mL
	172	Level 2	6 x 1 mL
	173	Level 3	6 x 1 mL
	172X	Trilevel MiniPak	3 x 1 mL

EXP 2019-05-31

LOT 38540

Level 1 38541
Level 2 38542
Level 3 38543

ENGLISH

INTENDED USE

Liquichek Diabetes Control is intended for use as an assayed quality control material to monitor the precision of laboratory testing procedures for the analytes listed in this package insert.

SUMMARY AND PRINCIPLE

The use of quality control materials is indicated as an objective assessment of the precision of methods and techniques in use and is an integral part of good laboratory practices. Three levels of control are available to allow performance monitoring within the clinical range. For customers in Germany, Quality control materials are required for assessment of laboratory performance as described in the "Guideline for Quality Assurance of Medical Laboratory Examinations following the German Medical Association" (RIL-BÄK regulation).

REAGENT

This product is prepared from human whole blood and contains preservatives and stabilizers. The control is provided in liquid form for convenience.

STORAGE AND STABILITY

This product will be stable until the expiration date when stored unopened at -10 to -70°C. The control may be stored unopened at 2 to 8°C for 6 months, but should not be used past the expiration date (note the date refrigerated temperature storage begins). Once the control is opened, all analytes will be stable for 14 days when stored tightly capped at 2 to 8°C.

Once thawed, do not refreeze the control; discard the remaining material.

This product is shipped under frozen conditions.

PROCEDURE

This product should be treated the same as patient specimens and run in accordance with instructions accompanying the instrument, kit or reagent being used.

Before sampling, allow the control to reach room temperature (18 to 25°C) and swirl gently to ensure homogeneity. Promptly replace the stopper and return to 2 to 8°C storage after each use.

Dispose of any discarded materials in accordance with the requirements of your local waste management authorities. In the event of damage to packaging, contact the local Bio-Rad Laboratories Sales Office or Bio-Rad Laboratories Technical Services.

LIMITATIONS

1. This product should not be used past the expiration date.
2. If there is evidence of microbial contamination or excessive turbidity in the product, discard the vial.
3. This product is not intended for use as a standard.
4. **Do not use this product if the packaging is damaged or the seal is broken.**

ASSIGNMENT OF VALUES

The mean values printed in this insert were derived from replicate analyses and are specific for this lot of product. The tests listed were performed by the manufacturer and/or independent laboratories using manufacturer supported reagents and a representative sampling of this lot of control. Individual laboratory means should fall within the corresponding acceptable range; however, laboratory means may vary from the listed values during the life of this control. Variations over time and between laboratories may be caused by differences in laboratory technique, instrumentation and reagents, or by manufacturer test method modifications. It is recommended that each laboratory establish its own means and acceptable ranges and use those provided only as guides.

Refer to www.qnet.com for insert update information.

SPECIFIC PERFORMANCE CHARACTERISTICS

This product is a stabilized liquid product manufactured under rigid quality control standards. To obtain consistent vial to vial assay values, the control requires proper storage and handling as described.

[I] Significant changes are highlighted.

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MEASUREMENT TARGET VALUES

The mean values given in this package insert were derived from replicate analyses and are specific for this product batch. The given tests were performed by the manufacturer and/or independent laboratories using reagents supported by the manufacturer and representative samples of this control batch. Individual laboratory mean values should fall within the acceptable range; laboratory mean values may vary from the given values over the life of the control. Variations over time and from laboratory to laboratory may be due to differences in laboratory technique, instruments and reagents, or to modifications of the manufacturer's test methods. We recommend that each laboratory establish its own mean values and acceptable ranges and only consider the given values as a guide.

Visit www.qcnet.com for current package insert information.

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This product is a stabilized liquid that is manufactured according to strict quality control standards. For reproducible analysis results from the vial to the vial, the control should be stored and handled correctly, as instructed.

[2] Important changes are highlighted!

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ENGLISH

- (1) NGSP units = %; IFCC unit = mmol/mol
- (2) The assigned values were obtained by using 1:300 dilution.
- ▲ Data is not available at this time. Please inquire.
- § The data required to establish the means and acceptable ranges for this assay were not obtained due to limited assignment participation. If your facility is interested in participating in the Value Assignment Program for this assay, please contact your local Bio-Rad office.
- ◆ INTERNATIONAL USE ONLY - The following section contains data for methods that are not available for diagnostic use in the United States.

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Diabetes/Hemoglobin Controls
Hematology Controls
Hepatitis B Response Controls
Immunology Controls
Immunology Protein Controls
Infectious Disease Controls
Specialty Infectious Disease Controls

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Abbott Laboratories, Abbott Park, Illinois

Arkray, Kyoto, Japan

Beckman Coulter, Inc., Brea, California

Beckman Coulter, Inc., Irving, Texas

Bio-Rad Laboratories, Inc., Hercules, California

Bio-Rad Laboratories, GmbH, Munich, Germany

HemoCue, Ängelholm, Sweden

Intopia, Gyeonggi-do, South Korea

Nihon Kohden Corporation, Tokyo, Japan

Ortho-Clinical Diagnostics, A Johnson and Johnson Company, Rochester, New York

Pointe Scientific, Inc., Canton, Michigan

PTS Diagnostics, Indianapolis, Indiana

Roche Diagnostics, GmbH, Penzberg, Germany

Roche Diagnostics, Indianapolis, Indiana

Sekisui Chemical Co. Ltd., Tokyo, Japan

Siemens Healthcare Diagnostics Inc., Newark, Delaware

Siemens Healthcare Diagnostics Inc., Tarrytown, New York

TOSOH Bioscience Inc., South San Francisco, California

TOSOH Bioscience Inc., Tokyo, Japan

Trinity Biotech, Kansas City, Missouri

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